

BAGE3 Antibody

Catalog No: #34476

Package Size: #34476-1 50ul #34476-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BAGE3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total BAGE3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human BAGE3.
Target Name	BAGE3
Other Names	Cancer/testis antigen 2.3; CT2.3;
Accession No.	Swiss-Prot: Q86Y29NCBI Gene ID: 85318
Uniprot	Q86Y29
GeneID	85318;
SDS-PAGE MW	12kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human testis tissue using BAGE3 antibody #34476.

Background

Unknown. Candidate gene encoding tumor antigens.

Ruault M., Eur. J. Hum. Genet. 10:833-840(2002).

Ruault M., Genomics 81:391-399(2003).

Note: This product is for in vitro research use only