

BAGE4 Antibody

Catalog No: #34477

Package Size: #34477-1 50ul #34477-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BAGE4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total BAGE4 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human BAGE4.
Target Name	BAGE4
Other Names	Cancer/testis antigen 2.4; CT2.4;
Accession No.	Swiss-Prot: Q86Y28NCBI Gene ID: 85317
Uniprot	Q86Y28
GeneID	85317;
SDS-PAGE MW	4kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using BAGE4 antibody #34477.

Background

Unknown. Candidate gene encoding tumor antigens.

Ruault M., Submitted (JAN-2001) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only