

EPB41L2 Antibody

Catalog No: #34478

Package Size: #34478-1 50ul #34478-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EPB41L2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total EPB41L2protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from Internal of human EPB41L2.
Target Name	EPB41L2
Other Names	4.1-G; 4.1G; Band 4.1-like protein 2; DKFZp781D1972; DKFZp781H1755
Accession No.	Swiss-Prot: O43491NCBI Gene ID: 2037
Uniprot	O43491
GeneID	2037;
SDS-PAGE MW	113kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, using EPB41L2 antibody #34478.

Background

2 isoforms of the human protein are produced by alternative splicing.

Parra M., Genomics 49:298-306(1998).

Liu J., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Mungall A.J., Nature 425:805-811(2003).

Note: This product is for in vitro research use only