

C9orf89 Antibody

Catalog No: #34479

Package Size: #34479-1 50ul #34479-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C9orf89 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total C9orf89 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human C9orf89.
Target Name	C9orf89
Other Names	Bcl10-interacting CARD protein; Bcl10-interacting protein with CARD; BINCA; BinCARD;
Accession No.	Swiss-Prot: Q96LW7NCBI Gene ID: 84270
Uniprot	Q96LW7
GeneID	84270;
SDS-PAGE MW	26kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human thyroid gland tissue using C9orf89 antibody #34479.

Background

Plays a role in inhibiting the effects of BCL10-induced activation of NF-kappa-B. May inhibit the phosphorylation of BCL10 in a CARD-dependent manner.

Humphray S.J., Nature 429:369-374(2004).

Woo H.-N., FEBS Lett. 578:239-244(2004).

Matsuoka S., Science 316:1160-1166(2007).

Note: This product is for in vitro research use only