

BNIP2 Antibody

Catalog No: #34480

Package Size: #34480-1 50ul #34480-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BNIP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total BNIP2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human BNIP2.
Target Name	BNIP2
Other Names	BCL2/adenovirus E1B 19-kDa protein-interacting 2; BCL2/adenovirus E1B 19kD interacting 2; BCL2/adenovirus E1B 19kD-interacting 2; BCL2/adenovirus E1B 19kDa interacting 2; BNIP-2
Accession No.	Swiss-Prot: Q12982NCBI Gene ID: 663
Uniprot	Q12982
GeneID	663;
SDS-PAGE MW	36kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, using BNIP2 antibody #34480.

Background

Implicated in the suppression of cell death. Interacts with the BCL-2 and adenovirus E1B 19 kDa proteins.

Boyd J.M., Cell 79:341-351(1994).

Rieder M.J., Submitted (APR-2003) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only