

GCNT6 Antibody

Catalog No: #34484

Package Size: #34484-1 50ul #34484-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GCNT6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total GCNT6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human GCNT6.
Target Name	GCNT6
Other Names	Beta-1; 3-galactosyl-O-glycosyl-glycoprotein beta-1; 6-N-acetylglucosaminyltransferase 6; EC 2.4.1.-;
Accession No.	Swiss-Prot: Q5T4J0NCBI Gene ID: 644378
Uniprot	Q5T4J0
GeneID	644378;
SDS-PAGE MW	45kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells and NIH/3T3 cells, using GCNT6 antibody #34484.

Background

Glycosyltransferase By similarity.

Mungall A.J., Nature 425:805-811(2003).

Note: This product is for in vitro research use only