

GLB1L3 Antibody

Catalog No: #34495

Package Size: #34495-1 50ul #34495-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GLB1L3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GLB1L3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human GLB1L3.
Target Name	GLB1L3
Other Names	Beta-galactosidase-1-like protein 3; EC 3.2.1.-;
Accession No.	Swiss-Prot: Q8NCI6NCBI Gene ID: 112937
Uniprot	Q8NCI6
GeneID	112937;
SDS-PAGE MW	70kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using GLB1L3 antibody #34495.

Background

Belongs to the glycosyl hydrolase 35 family.

Ota T., Nat. Genet. 36:40-45(2004).

Taylor T.D., Nature 440:497-500(2006).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only