

GUSBP1 Antibody

Catalog No: #34497

Package Size: #34497-1 50ul #34497-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GUSBP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GUSBP1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human GUSBP1.
Target Name	GUSBP1
Other Names	Beta-glucuronidase-like protein 1; SMA3-like protein;
Accession No.	Swiss-Prot: Q9NQQ0/Q15487/Q15486/Q15488/A6NI90/A8K1Q0/A8K3S8NCBI Gene ID: 153561
Uniprot	Q9NQQ0
GeneID	153561;
SDS-PAGE MW	18kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, using GUSBP1 antibody #34497.

Background

Sugano S., Nat. Genet. 36:40-45(2004).

Farmer A., Submitted (MAY-2003).

Rubin E.M., Nature 431:268-274(2004).

Note: This product is for in vitro research use only