

BAIAP2L1 Antibody

Catalog No: #34506

Package Size: #34506-1 50ul #34506-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BAIAP2L1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total BAIAP2L1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human BAIAP2L1.
Target Name	BAIAP2L1
Other Names	BAIAP2L1; BI2L1;
Accession No.	Swiss-Prot: Q9UHR4NCBI Gene ID: 55971
Uniprot	Q9UHR4
GeneID	55971;
SDS-PAGE MW	57kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HepG2 cells, using BAIAP2L1 antibody #34506.

Background

May function as adapter protein. Involved in the formation of clusters of actin bundles. Plays a role in the reorganization of the actin cytoskeleton in response to bacterial infection.

Hu R.-M., Proc. Natl. Acad. Sci. U.S.A. 97:9543-9548(2000).

Hillier L.W., Nature 424:157-164(2003).

Scherer S.W., Science 300:767-772(2003).

Note: This product is for in vitro research use only