

ABCC2 Antibody

Catalog No: #34526

Package Size: #34526-1 50ul #34526-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABCC2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total ABCC2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ABCC2.
Target Name	ABCC2
Other Names	Canalicular multispecific organic anion transporter 1; ATP-binding cassette sub-family C member 2; Canalicular multidrug resistance protein; Multidrug resistance-associated protein 2; ABCC2
Accession No.	Swiss-Prot: Q92887NCBI Gene ID: 1244
Uniprot	Q92887
GeneID	1244;
SDS-PAGE MW	174kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from RAW264.7 cells, using ABCC2 antibody #34526.

Background

Mediates hepatobiliary excretion of numerous organic anions. May function as a cellular cisplatin transporter.

Taniguchi K., Cancer Res. 56:4124-4129(1996).

Kool M., Submitted (FEB-1996) to the EMBL/GenBank/DDBJ databases.

Buechler M., J. Biol. Chem. 271:15091-15098(1996).

Note: This product is for in vitro research use only