

CLIP2 Antibody

Catalog No: #34530



Package Size: #34530-1 50ul #34530-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CLIP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CLIP2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human CLIP2.
Target Name	CLIP2
Other Names	CAP-Gly domain-containing linker protein 2; Cytoplasmic linker protein 115; CLIP-115; Cytoplasmic linker protein 2; Williams-Beuren syndrome chromosomal region 3 protein
Accession No.	Swiss-Prot: Q9UDT6NCBI Gene ID: 7461
Uniprot	Q9UDT6
GeneID	7461;
SDS-PAGE MW	120kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HepG2 cells and A549 cells, using CLIP2 antibody #34530.

Background

Seems to link microtubules to dendritic lamellar body (DLB), a membranous organelle predominantly present in bulbous dendritic appendages of neurons linked by dendrodendritic gap junctions. May operate in the control of brain-specific organelle translocations By similarity.

Ohara O., DNA Res. 4:53-59(1997).

Nakajima D., DNA Res. 9:99-106(2002).

Hillier L.W., Nature 424:157-164(2003).

Note: This product is for in vitro research use only