

CEP76 Antibody

Catalog No: #34569

Package Size: #34569-1 50ul #34569-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CEP76 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total CEP76 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human CEP76.
Target Name	CEP76
Other Names	C18orf9; centrosomal protein 76kDa; HsT1705;
Accession No.	Swiss-Prot: Q8TAP6NCBI Gene ID: 79959
Uniprot	Q8TAP6
GeneID	79959;
SDS-PAGE MW	74kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from 3T3 cells, using CEP76 antibody #34569.

Background

Centrosomal protein involved in regulation of centriole duplication. Required to limit centriole duplication to once per cell cycle by preventing centriole reduplication.

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Andersen J.S., Nature 426:570-574(2003)

Note: This product is for in vitro research use only