

CLIC3 Antibody

Catalog No: #34583

Package Size: #34583-1 50ul #34583-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLIC3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CLIC3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human CLIC3.
Target Name	CLIC3
Other Names	Chloride intracellular channel protein 3; CLIC3;
Accession No.	Swiss-Prot: O95833NCBI Gene ID: 9022
Uniprot	O95833
GeneID	9022;
SDS-PAGE MW	27kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, using CLIC3 antibody #34583.

Background

Can insert into membranes and form chloride ion channels. May participate in cellular growth control.

Humphray S.J., Nature 429:369-374(2004).

Mural R.J., Submitted (JUL-2005) to the EMBL/GenBank/DDBJ databases.

Qian Z., J. Biol. Chem. 274:1621-1627(1999).

Note: This product is for in vitro research use only