

## CDYLy2 Antibod

Catalog No: #34590

Package Size: #34590-1 50ul #34590-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | CDYLy2 Antibod                                                                                                                                     |
| Host Species          | Rabbit                                                                                                                                             |
| Clonality             | Polyclonal                                                                                                                                         |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                              |
| Applications          | IHC                                                                                                                                                |
| Species Reactivity    | Hu                                                                                                                                                 |
| Specificity           | The antibody detects endogenous levels of total CDYL2 protein.                                                                                     |
| Immunogen Type        | Peptide                                                                                                                                            |
| Immunogen Description | Synthesized peptide derived from N-terminal of humanCDYL2.                                                                                         |
| Target Name           | CDYLy2 Antibod                                                                                                                                     |
| Other Names           | Chromodomain Y-like protein 2; CDY-like 2;                                                                                                         |
| Accession No.         | Swiss-Prot: Q8N8U2NCBI Gene ID: 124359                                                                                                             |
| Uniprot               | Q8N8U2                                                                                                                                             |
| GeneID                | 124359;                                                                                                                                            |
| SDS-PAGE MW           | 57kd                                                                                                                                               |
| Concentration         | 1.0mg/ml                                                                                                                                           |
| Formulation           | Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at -20°C                                                                                                                                     |

## Application Details

Immunohistochemistry: 1:50~1:100

## Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using CDYL2 antibody #34590.

## Background

---

Ota T., Nat. Genet. 36:40-45(2004).

---

Note: This product is for in vitro research use only