

CSTF2T Antibody

Catalog No: #34599

Package Size: #34599-1 50ul #34599-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CSTF2T Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CSTF2T protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human CSTF2T.
Target Name	CSTF2T
Other Names	CF-1 64 kDa subunit; tau variant; cleavage stimulation factor 64 kDa subunit; tau variant; cleavage stimulation factor
Accession No.	Swiss-Prot: Q9H0L4NCBI Gene ID: 23283
Uniprot	Q9H0L4
GeneID	23283;
SDS-PAGE MW	64kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from MCF-7 cells and Jurkat cells, using CSTF2T antibody #34599.

Background

May play a significant role in AAUAAA-independent mRNA polyadenylation in germ cells. Directly involved in the binding to pre-mRNAs By similarity.

Dass B., Genomics 80:509-514(2002).

Wiemann S., Genome Res. 11:422-435(2001).

Beausoleil S.A., Proc. Natl. Acad. Sci. U.S.A. 101:12130-12135(2004).

Note: This product is for in vitro research use only