

CCDC102B Antibody

Catalog No: #34605

Package Size: #34605-1 50ul #34605-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	CCDC102B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CCDC102B protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human CCDC102B.
Target Name	CCDC102B
Other Names	Coiled-coil domain-containing protein 102B;
Accession No.	Swiss-Prot: Q68D86NCBI Gene ID: 79839
Uniprot	Q68D86
GeneID	79839;
SDS-PAGE MW	60kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from HuvEc cells and COLO cells, using CCDC102B antibody #34605.

Immunohistochemistry analysis of paraffin-embedded human placenta tissue, using CCDC102B antibody #34605.



Background

2 isoforms of the human protein are produced by alternative splicing.

Bechtel S., BMC Genomics 8:399-399(2007).

Nusbaum C., Nature 437:551-555(2005).

Ota T., Nat. Genet. 36:40-45(2004).

Note: This product is for in vitro research use only