

CCDC45 Antibody

Catalog No: #34606

Package Size: #34606-1 50ul #34606-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CCDC45 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CCDC45 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human CCDC45.
Target Name	CCDC45
Other Names	Coiled-coil domain-containing protein 45;
Accession No.	Swiss-Prot: Q96GE4NCBI Gene ID: 90799
Uniprot	Q96GE4
GeneID	90799;
SDS-PAGE MW	95kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HuvEc cells, using CCDC45 antibody #34606.

Background

Ota T., Nat. Genet. 36:40-45(2004).

Matsuoka S., Science 316:1160-1166(2007).

Dephoure N., Proc. Natl. Acad. Sci. U.S.A. 105:10762-10767(2008).

Note: This product is for in vitro research use only