

EEF1G Antibody

Catalog No: #34667



Package Size: #34667-1 50ul #34667-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	EEF1G Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total EEF1G protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human EEF1G.
Target Name	EEF1G
Other Names	eEF-1B gamma; EF-1-gamma; EF1G; elongation factor 1-gamma; eukaryotic translation elongation factor 1 gamma
Accession No.	Swiss-Prot: P26641NCBI Gene ID: 1937
Uniprot	P26641
GeneID	1937;
SDS-PAGE MW	55kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HepG2 cells and HUVEC cells, using EEF1G antibody #34667.

Background

Probably plays a role in anchoring the complex to other cellular components.

Sanders J., Nucleic Acids Res. 20:5907-5910(1992).

Kumabe T., Nucleic Acids Res. 20:2598-2598(1992).

Gevaert K., Nat. Biotechnol. 21:566-569(2003).

Note: This product is for in vitro research use only