

FKBPL Antibody

Catalog No: #34682

Package Size: #34682-1 50ul #34682-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FKBPL Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total FKBPL protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human FKBPL.
Target Name	FKBPL
Other Names	DIR1; FK506 binding protein like; WAF-1/CIP1 stabilizing 39; WISP39;
Accession No.	Swiss-Prot: Q9UIM3NCBI Gene ID: 63943
Uniprot	Q9UIM3
GeneID	63943;
SDS-PAGE MW	38kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, using FKBPL antibody #34682.

Background

May be involved in response to X-ray. Regulates p21 protein stability by binding to Hsp90 and p21.

Robson T., Radiat. Res. 152:451-461(1999).

Ota T., Nat. Genet. 36:40-45(2004).

Xie T., Genome Res. 13:2621-2636(2003).

Note: This product is for in vitro research use only