

FOXB2 Antibody

Catalog No: #34684

Package Size: #34684-1 50ul #34684-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FOXB2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total FOXB2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human FOXB2.
Target Name	FOXB2
Other Names	forkhead box B2;
Accession No.	Swiss-Prot: Q5VYV0NCBI Gene ID: 442425
Uniprot	Q5VYV0
GeneID	442425;
SDS-PAGE MW	46kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HeLa cells, LOVO cells and A549 cells, using FOXB2 antibody #34684.

Background

Transcription factor Probable.

Humphray S.J., Nature 429:369-374(2004).

Note: This product is for in vitro research use only