

FYCO1 Antibody

Catalog No: #34690

Package Size: #34690-1 50ul #34690-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FYCO1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FYCO1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from Internal of human FYCO1.
Target Name	FYCO1
Other Names	DKFZp779K1152; FLJ13335; FYCO1; FYVE and coiled-coil domain containing 1; FYVE and coiled-coil domain-containing protein 1
Accession No.	Swiss-Prot: Q9BQS8NCBI Gene ID: 79443
Uniprot	Q9BQS8
GeneID	79443;
SDS-PAGE MW	168kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from JurKat cells, using FYCO1 antibody #34690.

Background

May mediate microtubule plus end-directed vesicle transport.

Kiss H., Eur. J. Hum. Genet. 10:52-61(2002).

Bechtel S., BMC Genomics 8:399-399(2007).

Muzny D.M., Nature 440:1194-1198(2006).

Note: This product is for in vitro research use only