

GPRIN1 Antibody

Catalog No: #34691

Package Size: #34691-1 50ul #34691-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPRIN1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GPRIN1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human GPRIN1.
Target Name	GPRIN1
Other Names	G protein-regulated inducer of neurite outgrowth 1; KIAA1893;
Accession No.	Swiss-Prot: Q7Z2K8NCBI Gene ID: 114787
Uniprot	Q7Z2K8
GeneID	114787;
SDS-PAGE MW	102kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HT-29 cells, using GPRIN1 antibody #34691.

Background

May be involved in neurite outgrowth By similarity.

Nagase T., DNA Res. 8:179-187(2001).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

The German cDNA consortium; Submitted (JUL-2002) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only