

TUBGCP6 Antibody

Catalog No: #34698

Package Size: #34698-1 50ul #34698-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	TUBGCP6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TUBGCP6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human TUBGCP6.
Target Name	TUBGCP6
Other Names	gamma-tubulin complex component 6; GCP-6; GCP6; KIAA1669;
Accession No.	Swiss-Prot: Q96RT7NCBI Gene ID: 85378
Uniprot	Q96RT7
GeneID	85378;
SDS-PAGE MW	200kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details


Western blotting: 1:500~1:3000

Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from COLO cells, using TUBGCP6 antibody #34698.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using TUBGCP6 antibody #34698.

Background

Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome.

Murphy S.M., Mol. Biol. Cell 12:3340-3352(2001).

Dunham I., Nature 402:489-495(1999).

Hirosawa M., DNA Res. 8:1-9(2001).

Note: This product is for in vitro research use only