

GNL1 Antibody

Catalog No: #34731



Package Size: #34731-1 50ul #34731-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GNL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total GNL1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human GNL1.
Target Name	GNL1
Other Names	GTP-binding protein HSR1; guanine nucleotide binding protein-like 1; guanine nucleotide-binding protein-like 1; HSR1; HSR1 GTP-binding protein
Accession No.	Swiss-Prot: P36915NCBI Gene ID: 2794
Uniprot	P36915
GeneID	2794;
SDS-PAGE MW	69kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, LOVO cells and K562 cells, using GNL1 antibody #34731.

Background

Possible regulatory or functional link with the histocompatibility cluster.

Vernet C., Mamm. Genome 5:100-105(1994).

Kalnina N., Submitted (MAY-2003) to the EMBL/GenBank/DDBJ databases.

Mungall A.J., Nature 425:805-811(2003).

Note: This product is for in vitro research use only