

HNRPLL Antibody

Catalog No: #34733

Package Size: #34733-1 50ul #34733-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HNRPLL Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total HNRPLL protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human HNRPLL.
Target Name	HNRPLL
Other Names	BLOCK24 protein; heterogeneous nuclear ribonucleoprotein L-like; HNRLL; SRRF; stromal RNA-regulating factor
Accession No.	Swiss-Prot: Q8WVV9NCBI Gene ID: 92906
Uniprot	Q8WVV9
GeneID	92906;
SDS-PAGE MW	60kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from COLO cells and HUVEC cells, using HNRPLL antibody #34733.

Background

RNA-binding protein that functions as regulator of alternative splicing for multiple target mRNAs, including PTPRC/CD45 and STAT5A. Required for alternative splicing of PTPRC.

Gorry M.C., Submitted (DEC-2001) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

The German cDNA consortium; Submitted (SEP-2004) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only