

HLAH Antibody

Catalog No: #34734

Package Size: #34734-1 50ul #34734-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HLAH Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total HLAH protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human HLAH.
Target Name	HLAH
Other Names	HLA class I histocompatibility antigen; alpha chain H precursor; HLA-12.4; HLA-AR; HLAH
Accession No.	Swiss-Prot: P01893NCBI Gene ID: 3136
Uniprot	P01893
GeneID	3136;
SDS-PAGE MW	40kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from LOVO cells, using HLAH antibody #34734.

Background

Involved in the presentation of foreign antigens to the immune system.

Malissen M., Proc. Natl. Acad. Sci. U.S.A. 79:893-897(1982).

Zahedi R.P., J. Proteome Res. 7:526-534(2008).

Note: This product is for in vitro research use only