

HOXA6 Antibody

Catalog No: #34741

Package Size: #34741-1 50ul #34741-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HOXA6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total HOXA6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human HOXA6.
Target Name	HOXA6
Other Names	Homeobox protein Hox-A6; Hox-1B; HOXA6; HOX1B;
Accession No.	Swiss-Prot: P31267NCBI Gene ID: 3203
Uniprot	P31267
GeneID	3203;
SDS-PAGE MW	26kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HT-29 cells and K562 cells, using HOXA6 antibody #34741.

Background

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Hillier L.W., Nature 424:157-164(2003).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Boncinelli E., Genome 31:745-756(1989).

Note: This product is for in vitro research use only