

USP50 Antibody

Catalog No: #34748

Package Size: #34748-1 50ul #34748-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	USP50 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total USP50 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human USP50.
Target Name	USP50
Other Names	Inactive ubiquitin carboxyl-terminal hydrolase 50; Inactive ubiquitin-specific peptidase 50;
Accession No.	Swiss-Prot: Q70EL3NCBI Gene ID: 373509
Uniprot	Q70EL3
GeneID	373509;
SDS-PAGE MW	38kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from A549 cells, using USP50 antibody #34748.

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using USP50 antibody #34748.



Background

Has no peptidase activity.

Quesada V., m. Biophys. Res. Commun. 314:54-62(2004).

Puente X.S., Genome Res. 14:609-622(2004).

Note: This product is for in vitro research use only