

IQCB1 Antibody

Catalog No: #34750

Package Size: #34750-1 50ul #34750-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IQCB1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total IQCB1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human IQCB1.
Target Name	IQCB1
Other Names	IQ motif containing B1; NPHP5; PIQ; SLSN5;
Accession No.	Swiss-Prot: Q15051NCBI Gene ID: 9657
Uniprot	Q15051
GeneID	9657;
SDS-PAGE MW	69kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, using IQCB1 antibody #34750.

Background

Involved in cillogenesis By similarity.

Otto E.A., Nat. Genet. 37:282-288(2005).

Nomura N., DNA Res. 1:27-35(1994).

Shichijo S., Submitted (MAY-2001) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only