

KLHL29 Antibody

Catalog No: #34751

Package Size: #34751-1 50ul #34751-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KLHL29 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total KLHL29 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human KLHL29.
Target Name	KLHL29
Other Names	Kelch-like protein 29; Kelch repeat and BTB domain-containing protein 9;
Accession No.	Swiss-Prot: Q96CT2NCBI Gene ID: 114818
Uniprot	Q96CT2
GeneID	114818;
SDS-PAGE MW	72kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue using KLHL29 antibody #34751.

Background

Bechtel S., BMC Genomics 8:399-399(2007).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Nagase T., DNA Res. 8:179-187(2001).

Note: This product is for in vitro research use only