

LILRA1 Antibody

Catalog No: #34757

Package Size: #34757-1 50ul #34757-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	LILRA1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total LILRA1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human LILRA1.
Target Name	LILRA1
Other Names	CD85l; leukocyte immunoglobulin-like receptor 6; leukocyte immunoglobulin-like receptor subfamily A member 1; LIR-6; LIR6
Accession No.	Swiss-Prot: O75019NCBI Gene ID: 11024
Uniprot	O75019
GeneID	11024;
SDS-PAGE MW	53kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunofluorescence: 1:100~1:500

Images



Western blot analysis of extracts from HeLa cells and MCF-7 cells, using LILRA1 antibody #34757.

Immunofluorescence analysis of A549 cells, using LILRA1 antibody #34757.



Background

May act as receptor for class I MHC antigens.

Borges L., J. Immunol. 159:5192-5196(1997).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only