

LILRA2 Antibody

Catalog No: #34758

Package Size: #34758-1 50ul #34758-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LILRA2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total LILRA2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human LILRA2.
Target Name	LILRA2
Other Names	Leukocyte immunoglobulin-like receptor subfamily A member 2 [Precursor]; Leukocyte immunoglobulin-like receptor 7; LIR-7; Immunoglobulin-like transcript 1; ILT-1
Accession No.	Swiss-Prot: Q8N149NCBI Gene ID: 11027
Uniprot	Q8N149
GeneID	11027;
SDS-PAGE MW	53kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HepG2 cells, using LILRA2 antibody #34758.

Background

Does not bind class I MHC antigens.

Borges L., J. Immunol. 159:5192-5196(1997).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only