

LDLRAD2 Antibody

Catalog No: #34764

Package Size: #34764-1 50ul #34764-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	LDLRAD2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total LDLRAD2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human LDLRAD2.
Target Name	LDLRAD2
Other Names	Idlrad2; low density lipoprotein receptor class a domain containing 2; lrad2;
Accession No.	Swiss-Prot: Q5SZI1NCBI Gene ID: 401944
Uniprot	Q5SZI1
GeneID	401944;
SDS-PAGE MW	29kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, LOVO cells and 3T3 using LDLRAD2 antibody #34764.

Background

Belongs to the LDLR family.

Ota T., Nat. Genet. 36:40-45(2004).

Gregory S.G., Nature 441:315-321(2006).

Mural R.J., Submitted (JUL-2005) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only