

LRP10 Antibody

Catalog No: #34766

Package Size: #34766-1 50ul #34766-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LRP10 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total LRP10 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human LRP10.
Target Name	LRP10
Other Names	Low-density lipoprotein receptor-related protein 10;
Accession No.	Swiss-Prot: Q7Z4F1NCBI Gene ID: 26020
Uniprot	Q7Z4F1
GeneID	26020;
SDS-PAGE MW	76kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue using LRP10 antibody #34766.

Background

Probable receptor, which is involved in the internalization of lipophilic molecules and/or signal transduction. May be involved in the uptake of lipoprotein APOE in liver By similarity.

Clark H.F., Genome Res. 13:2265-2270(2003).

Li W.B., Submitted (JAN-2003) to the EMBL/GenBank/DDBJ databases.

Bechtel S., BMC Genomics 8:399-399(2007).

Note: This product is for in vitro research use only