

LAMP3 Antibody

Catalog No: #34772

Package Size: #34772-1 50ul #34772-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LAMP3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total LAMP3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human LAMP3.
Target Name	LAMP3
Other Names	Lysosome-associated membrane glycoprotein 3 ; LAMP-3; Lysosomal-associated membrane protein 3; DC-lysosome-associated membrane glycoprotein; DC LAMP
Accession No.	Swiss-Prot: Q9UQV4NCBI Gene ID: 27074
Uniprot	Q9UQV4
GeneID	27074;
SDS-PAGE MW	44kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using LAMP3 antibody #34772.

Background

May play a role in dendritic cell function and in adaptive immunity.

Ozaki K., Cancer Res. 58:3499-3503(1998).

Banchereau J., Immunity 9:325-336(1998).

Note: This product is for in vitro research use only