

MIA2 Antibody

Catalog No: #34780

Package Size: #34780-1 50ul #34780-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MIA2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total MIA2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human MIA2.
Target Name	MIA2
Other Names	Melanoma inhibitory activity protein 2;
Accession No.	Swiss-Prot: Q96PC5NCBI Gene ID: 117153
Uniprot	Q96PC5
GeneID	4253;
SDS-PAGE MW	62kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HuvEc cells, using MIA2 antibody #34780.

Background

May play a role in the pathophysiology of liver disease and may serve as a marker of liver damage.

Bosserhoff A.K., J. Biol. Chem. 278:15225-15231(2003).

Heilig R., Nature 421:601-607(2003).

Mural R.J., Submitted (JUL-2005) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only