

MIA Antibody

Catalog No: #34781

Package Size: #34781-1 50ul #34781-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MIA Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total MIA protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human MIA.
Target Name	MIA
Other Names	Melanoma-derived growth regulatory protein; Melanoma inhibitory activity;
Accession No.	Swiss-Prot: Q16674NCBI Gene ID: 8190
Uniprot	Q16674
GeneID	8190;
SDS-PAGE MW	15kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human heart tissue using MIA antibody #34781.

Background

Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.

Blesch A., Cancer Res. 54:5695-5701(1994).

Bosserhoff.,J. Biol. Chem. 271:490-495(1996).

Note: This product is for in vitro research use only