

GSDMC Antibody

Catalog No: #34782

Package Size: #34782-1 50ul #34782-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GSDMC Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GSDMC protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human GSDMC.
Target Name	GSDMC
Other Names	Melanoma-derived leucine zipper-containing extranuclear factor;
Accession No.	Swiss-Prot: Q9BYG8NCBI Gene ID: 56169
Uniprot	Q9BYG8
GeneID	56169;
SDS-PAGE MW	58kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from 293 cell, using GSDMC antibody #34782.

Background

Belongs to the gasdermin family.

Watabe K., Jpn. J. Cancer Res. 92:140-151(2001).

Nusbaum C., Nature 439:331-335(2006).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only