

NSG2 Antibody

Catalog No: #34861

Package Size: #34861-1 50ul #34861-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NSG2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total NSG2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human NSG2.
Target Name	NSG2
Other Names	Neuron-specific protein family member 2; Protein p19; HMP19;
Accession No.	Swiss-Prot: Q9Y328NCBI Gene ID: 51617
Uniprot	Q9Y328
GeneID	51617;
SDS-PAGE MW	19kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue using NSG2 antibody #34861.

Background

Belongs to the NSG family.

Jin J., Proc. Natl. Acad. Sci. U.S.A. 97:9543-9548(2000).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only