

NCOA7 Antibody

Catalog No: #34871

Package Size: #34871-1 50ul #34871-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NCOA7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total NCOA7 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human NCOA7.
Target Name	NCOA7
Other Names	ERAP140; ESNA1; estrogen nuclear receptor coactivator 1; estrogen receptor associated protein 140 kDa; nuclear receptor coactivator 7
Accession No.	Swiss-Prot: Q8NI08NCBI Gene ID: 135112
Uniprot	Q8NI08
GeneID	135112;
SDS-PAGE MW	106kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from LOVO cells, 3T3 cells and RAW264.7 cells, using NCOA7 antibody #34871.

Background

Enhances the transcriptional activities of several nuclear receptors. Involved in the coactivation of different nuclear receptors, such as ESR1, THRB, PPARG and RARA.

Shao W., Mol. Cell. Biol. 22:3358-3372(2002).

Ohira M., Cancer Lett. 197:63-68(2003).

Ota T., Nat. Genet. 36:40-45(2004).

Note: This product is for in vitro research use only