

NXF3 Antibody

Catalog No: #34872

Package Size: #34872-1 50ul #34872-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NXF3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NXF3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human NXF3.
Target Name	NXF3
Other Names	nuclear RNA export factor 3; TAP-like 3; TAPL-3; TAPL3;
Accession No.	Swiss-Prot: Q9H4D5NCBI Gene ID: 56000
Uniprot	Q9H4D5
GeneID	56000;
SDS-PAGE MW	60kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from MCF-7 cells, using NXF3 antibody #34872.

Background

May function as a tissue-specific nuclear mRNA export factor.

Herold A., Mol. Cell. Biol. 20:8996-9008(2000).

Yang J., Mol. Cell 8:397-406(2001).

The MGC Project Team; Genome Res. 14:2121-2127(2004)

Note: This product is for in vitro research use only