

NUSAP1 Antibody

Catalog No: #34873

Package Size: #34873-1 50ul #34873-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NUSAP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NUSAP1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human NUSAP1.
Target Name	NUSAP1
Other Names	ANKT; BM037; LNP; nucleolar and spindle associated protein 1; nucleolar protein ANKT
Accession No.	Swiss-Prot: Q9BXS6NCBI Gene ID: 51203
Uniprot	Q9BXS6
GeneID	51203;
SDS-PAGE MW	49kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HT-29 cells, using NUSAP1 antibody #34873.

Background

Microtubule-associated protein with the capacity to bundle and stabilize microtubules By similarity. May associate with chromosomes and promote the organization of mitotic spindle microtubules around them.

Qu X., Submitted (JUL-2000) to the EMBL/GenBank/DDBJ databases.

Sato H., Submitted (OCT-2000) to the EMBL/GenBank/DDBJ databases.

Yu Y., Genome Res. 11:1392-1403(2001)

Note: This product is for in vitro research use only