

NOC3L Antibody

Catalog No: #34874

Package Size: #34874-1 50ul #34874-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NOC3L Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total NOC3L protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human NOC3L.
Target Name	NOC3L
Other Names	AD24; C10orf117; FAD24; FLJ12820; nucleolar complex associated 3 homolog
Accession No.	Swiss-Prot: Q8WTT2NCBI Gene ID: 64318
Uniprot	Q8WTT2
GeneID	64318;
SDS-PAGE MW	93kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, using NOC3L antibody #34874.

Background

May be required for adipogenesis By similarity.

Tominaga K., J. Cell Sci. 117:6217-6226(2004).

Deloukas P., Nature 429:375-381(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004)

Note: This product is for in vitro research use only