

URB1 Antibody

Catalog No: #34877

Package Size: #34877-1 50ul #34877-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	URB1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total URB1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human URB1.
Target Name	URB1
Other Names	Nucleolar pre-ribosomal-associated protein 1; Protein URB1;
Accession No.	Swiss-Prot: O60287NCBI Gene ID: 9875
Uniprot	O60287
GeneID	9875;
SDS-PAGE MW	255kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue using URB1 antibody #34877.

Background

The German cDNA consortium; Submitted (JUN-2005) to the EMBL/GenBank/DDBJ databases.

Slavov D., Gene 247:215-232(2000).

Nagase T., DNA Res. 5:31-39(1998).

Note: This product is for in vitro research use only