

NR4A3 Antibody

Catalog No: #34884

Package Size: #34884-1 50ul #34884-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NR4A3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NR4A3 protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human NR4A3.
Target Name	NR4A3
Other Names	Nuclear receptor subfamily 4 group A member 3; Nuclear hormone receptor NOR-1; Neuron-derived orphan receptor 1; Mitogen-induced nuclear orphan receptor;
Accession No.	Swiss-Prot: Q92570NCBI Gene ID: 8013
Uniprot	Q92570
GeneID	8013;
SDS-PAGE MW	65kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from COLO cells, using NR4A3 antibody #34884.

Background

Binds to the B1A response-element.

Ohkura N., Biochim. Biophys. Acta 1308:205-214(1996).

Hedvat C.V., Mol. Endocrinol. 9:1692-1700(1995).

Clark J., Oncogene 12:229-235(1996).

Note: This product is for in vitro research use only