

OXR1 Antibody

Catalog No: #34885

Package Size: #34885-1 50ul #34885-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OXR1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total OXR1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human OXR1.
Target Name	OXR1
Other Names	Oxidation resistance protein 1;
Accession No.	Swiss-Prot: Q8N573NCBI Gene ID: 55074
Uniprot	Q8N573
GeneID	55074;
SDS-PAGE MW	85kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue using OXR1 antibody #34885.

Background

May be involved in protection from oxidative damage.

Ohira M., Cancer Lett. 197:63-68(2003).

The German cDNA consortium; Submitted (JUN-2005) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only