

PHF3 Antibody

Catalog No: #34899

Package Size: #34899-1 50ul #34899-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	PHF3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total PHF3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human PHF3.
Target Name	PHF3
Other Names	PHD finger protein 3;
Accession No.	Swiss-Prot: Q92576NCBI Gene ID: 23469
Uniprot	Q92576
GeneID	23469;
SDS-PAGE MW	230kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000


Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from HeLa cells, using PHF3 antibody #34899.

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using PHF3 antibody #34899.



Background

This gene encodes a member of a PHD finger-containing gene family. This gene may function as a transcription factor and may be involved in glioblastomas development. Alternative splicing results in multiple transcript variants.

Fischer U., Cytogenet. Cell Genet. 94:131-136(2001).

Nagase T., DNA Res. 3:321-329(1996).

Nakajima D., DNA Res. 9:99-106(2002).

Note: This product is for in vitro research use only