

N4BP2L2 Antibody

Catalog No: #34906

Package Size: #34906-1 50ul #34906-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	N4BP2L2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total N4BP2L2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human N4BP2L2.
Target Name	N4BP2L2
Other Names	NEDD4-binding protein 2-like 2; Phosphonoformate immuno-associated protein 5; N4BP2L2; CG005; PFAAP5
Accession No.	Swiss-Prot: Q92802NCBI Gene ID: 10443
Uniprot	Q92802
GeneID	10443;
SDS-PAGE MW	67kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from LOVO cells, using N4BP2L2 antibody #34906.

Background

3 isoforms of the human protein are produced by alternative splicing.

Couch F.J., Genomics 36:86-99(1996).

Liu Y., Submitted (JUL-2002) to the EMBL/GenBank/DDBJ databases.

Kalnina N., Submitted (MAY-2003) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only