

PRPF39 Antibody

Catalog No: #34924

Package Size: #34924-1 50ul #34924-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PRPF39 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total PRPF39 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human PRPF39.
Target Name	PRPF39
Other Names	flj11128; flj20666; flj45460; mgc149842; mgc149843
Accession No.	Swiss-Prot: Q86UA1NCBI Gene ID: 55015
Uniprot	Q86UA1
GeneID	55015;
SDS-PAGE MW	74kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, using PRPF39 antibody #34924.

Background

Involved in pre-mRNA splicing By similarity.

Ota T.,Nat. Genet. 36:40-45(2004).

The MGC Project Team;Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only