

RABEP2 Antibody

Catalog No: #34947

Package Size: #34947-1 50ul #34947-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RABEP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total RABEP2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RABEP2.
Target Name	RABEP2
Other Names	rab GTPase-binding effector protein 2; rabaptin-5beta; rabe2; rabpt5b;
Accession No.	Swiss-Prot: Q9H5N1NCBI Gene ID: 79874
Uniprot	Q9H5N1
GeneID	79874;
SDS-PAGE MW	64kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HeLa cells, Jurkat cells, 293 cells and HUVEC cells, using RABEP2 antibody #34947.

Background

Plays a role in membrane trafficking and in homotypic early endosome fusion.

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Gournier H., EMBO J. 17:1930-1940(1998).

Note: This product is for in vitro research use only